

Human BLNK/B-cell linker protein ELISA Kit

Catalog No: #EK19769

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BLNK/B-cell linker protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human B-cell linker protein
Target Name	B-cell linker protein
Other Names	B-cell adapter containing a SH2 domain protein, B-cell adapter containing a Src homology 2 domain protein, Cytoplasmic adapter protein, Src homology 2 domain-containing leukocyte protein of 65 kDa, SLP-65, BASH, SLP65
Accession No.	Q8WV28
Uniprot	Q8WV28
GeneID	29760;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.084ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BLNK

Note: This product is for in vitro research use only