

Mouse Blnk/B-cell linker protein ELISA Kit

Catalog No: #EK19770

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Blnk/B-cell linker protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse B-cell linker protein
Target Name	B-cell linker protein
Other Names	B-cell adapter containing a SH2 domain protein,B-cell adapter containing a Src homology 2 domain protein,Cytoplasmic adapter protein,Lymphocyte antigen 57,Src homology 2 domain-containing leukocyte protein of 65 kDa,Slp-65,Bash,Ly57,Slp65
Accession No.	Q9QUN3
Uniprot	Q9QUN3
GeneID	17060;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:41.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Blnk

Note: This product is for in vitro research use only