

## Rat Blnk/B-cell linker protein ELISA Kit

Catalog No: #EK19771

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Blnk/B-cell linker protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat B-cell linker protein
Target Name	B-cell linker protein
Other Names	Cytoplasmic adapter protein
Accession No.	Q4KM52
Uniprot	Q4KM52
GeneID	499356;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.088ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Blnk

Note: This product is for in vitro research use only