

## RIOK2 Antibody

Catalog No: #44685



Package Size: #44685-1 50ul #44685-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RIOK2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	RIOK2antibody detects endogenous levels of total?RIOK2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RIOK2
Target Name	RIOK2
Other Names	FLJ11159 antibody OTTHUMP00000158911 antibody OTTHUMP00000222868 antibody RIO kinase 2 (yeast) antibody RIO kinase 2 antibody RIO2 antibody RIOK2 antibody RIOK2_HUMAN antibody Serine/threonine protein kinase RIO2 antibody Serine/threonine-protein kinase rio2 antibody
Accession No.	Swiss-Prot#: Q9BVS4NCBI Gene ID:55781
Uniprot	Q9BVS4
GeneID	55781;
Calculated MW	63 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only