

Human ANK2/Ankyrin-2 ELISA Kit

Catalog No: #EK19773

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ANK2/Ankyrin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Ankyrin-2
Target Name	Ankyrin-2
Other Names	ANK-2, Ankyrin-B, Brain ankyrin, Non-erythroid ankyrin
Accession No.	Q01484
Uniprot	Q01484
GeneID	287;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.41ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ANK2

Note: This product is for in vitro research use only