

## Mouse Ank2/Ankyrin-2 ELISA Kit

Catalog No: #EK19774

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ank2/Ankyrin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Ankyrin-2
Target Name	Ankyrin-2
Other Names	ANK-2,Ankyrin-B,Brain ankyrin
Accession No.	Q8C8R3
Uniprot	Q8C8R3
GeneID	109676;

## Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.411ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Ank2

Note: This product is for in vitro research use only