

## Mouse Acta1/Actin, alpha skeletal muscle ELISA Kit

Catalog No: #EK19776

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Acta1/Actin, alpha skeletal muscle ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Actin, alpha skeletal muscle
Target Name	Actin, alpha skeletal muscle
Other Names	Alpha-actin-1,Acta
Accession No.	P68134
Uniprot	P68134
GeneID	11459;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.092ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Acta1

Note: This product is for in vitro research use only