

Human CLCF1/Cardiotrophin-like cytokine factor 1 ELISA Kit



Catalog No: #EK19778

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CLCF1/Cardiotrophin-like cytokine factor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cardiotrophin-like cytokine factor 1
Target Name	Cardiotrophin-like cytokine factor 1
Other Names	B-cell-stimulating factor 3,BSF-3,Novel neurotrophin-1,NNT-1,BSF3,CLC,NNT1
Accession No.	Q9UBD9
Uniprot	Q9UBD9
GeneID	23529;

Application Details

Assay Type: Sandwich
Detect Range:15.6-1000pg/mL
Sensitivity:7.8pg/mL
Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul
Assay Time:3-5h
Detection wavelength:450 nm
gene_full_name:CLCF1

Note: This product is for in vitro research use only