

Mouse Clcf1/Cardiotrophin-like cytokine factor 1 ELISA Kit

Catalog No: #EK19779

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Clcf1/Cardiotrophin-like cytokine factor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cardiotrophin-like cytokine factor 1
Target Name	Cardiotrophin-like cytokine factor 1
Other Names	B-cell-stimulating factor 3,BSF-3,Novel neurotrophin-1,NNT-1,Bsf3,Nnt1
Accession No.	Q9QZM3
Uniprot	Q9QZM3
GeneID	56708;

Application Details

Assay Type: SandwichDetect Range:15.6-1000pg/mL Sensitivity:7.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Clcf1

Note: This product is for in vitro research use only