

Human CCR7/C-C chemokine receptor type 7 ELISA Kit

Catalog No: #EK19780

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CCR7/C-C chemokine receptor type 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C chemokine receptor type 7
Target Name	C-C chemokine receptor type 7
Other Names	C-C CKR-7,BLR2,CDw197,Epstein-Barr virus-induced G-protein coupled receptor 1,EBI1,MIP-3 beta receptor,CD197,CMKBR7,EBI1,EVI1
Accession No.	P32248
Uniprot	P32248
GeneID	1236;

Application Details

Assay Type: SandwichDetect Range:31.2-2000pg/mL Sensitivity:15.78pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CCR7

Note: This product is for in vitro research use only