

Mouse Ccr7/C-C chemokine receptor type 7 ELISA Kit

Catalog No: #EK19781

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ccr7/C-C chemokine receptor type 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C chemokine receptor type 7
Target Name	C-C chemokine receptor type 7
Other Names	C-C CKR-7, Epstein-Barr virus-induced G-protein coupled receptor 1, EBI1, MIP-3 beta receptor, CD197, Cmkbr7, Ebi1, Ebi1h
Accession No.	P47774
Uniprot	P47774
GeneID	12775;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.94ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ccr7

Note: This product is for in vitro research use only