

## Human CCL27/C-C motif chemokine 27 ELISA Kit

Catalog No: #EK19782

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CCL27/C-C motif chemokine 27 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C motif chemokine 27
Target Name	C-C motif chemokine 27
Other Names	CC chemokine ILC,Cutaneous T-cell-attracting chemokine,CTACK,ESkine,IL-11 R-alpha-locus chemokine,Skinkine,Small-inducible cytokine A27,ILC,SCYA27
Accession No.	Q9Y4X3
Uniprot	Q9Y4X3
GeneID	10850;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/ml Sensitivity:42.5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CCL27

Note: This product is for in vitro research use only