

Mouse Ccl27/C-C motif chemokine 27 ELISA Kit

Catalog No: #EK19783

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ccl27/C-C motif chemokine 27 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C motif chemokine 27
Target Name	C-C motif chemokine 27
Other Names	CC chemokine ILC,Cutaneous T-cell-attracting chemokine,CTACK,ESkine,IL-11 R-alpha-locus chemokine,ALP,Skinkine,Small-inducible cytokine A27,Ilc,Scya27
Accession No.	Q9Z1X0
Uniprot	Q9Z1X0
GeneID	100039863;100039939;100040048;100861978;20301;

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.31ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ccl27

Note: This product is for in vitro research use only