

Human AGTR1/Type-1 angiotensin II receptor ELISA Kit



Catalog No: #EK19784

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Human AGTR1/Type-1 angiotensin II receptor ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Hu |
| Specificity | Natural and recombinant Human Type-1 angiotensin II receptor |
| Target Name | Type-1 angiotensin II receptor |
| Other Names | AT1AR,AT1BR,Angiotensin II type-1 receptor,AT1,AGTR1A,AGTR1B,AT2R1,AT2R1B |
| Accession No. | P30556 |
| Uniprot | P30556 |
| GeneID | 185; |

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.3ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AGTR1

Note: This product is for in vitro research use only