

Mouse Agtr1a/Type-1A angiotensin II receptor ELISA Kit

Catalog No: #EK19785

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Agtr1a/Type-1A angiotensin II receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Type-1A angiotensin II receptor
Target Name	Type-1A angiotensin II receptor
Other Names	Angiotensin II type-1A receptor,AT1A,Agtr1
Accession No.	P29754
Uniprot	P29754
GeneID	11607;

Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.413ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Agtr1a

Note: This product is for in vitro research use only