

## Chicken CS/Citrate synthase, mitochondrial ELISA Kit

Catalog No: #EK19786

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken CS/Citrate synthase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Citrate synthase, mitochondrial
Target Name	Citrate synthase, mitochondrial
Other Names	Citrate (Si)-synthase
Accession No.	P23007
Uniprot	P23007

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CS

Note: This product is for in vitro research use only