

Human CS/Citrate synthase, mitochondrial ELISA Kit

Catalog No: #EK19787

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CS/Citrate synthase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Citrate synthase, mitochondrial
Target Name	Citrate synthase, mitochondrial
Other Names	Citrate (Si)-synthase
Accession No.	O75390
Uniprot	O75390
GeneID	1431;

Application Details

Assay Type: SandwichDetect Range:31.2-2000pg/mL Sensitivity:15.79pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CS

Note: This product is for in vitro research use only