

Mouse Cs/Citrate synthase, mitochondrial ELISA Kit

Catalog No: #EK19788

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cs/Citrate synthase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Citrate synthase, mitochondrial
Target Name	Citrate synthase, mitochondrial
Other Names	Citrate (Si)-synthase
Accession No.	Q9CZU6
Uniprot	Q9CZU6
GeneID	12974;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cs

Note: This product is for in vitro research use only