

Rat Cs/Citrate synthase, mitochondrial ELISA Kit

Catalog No: #EK19789

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cs/Citrate synthase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Citrate synthase, mitochondrial
Target Name	Citrate synthase, mitochondrial
Other Names	Citrate (Si)-synthase
Accession No.	Q8VHF5
Uniprot	Q8VHF5
GeneID	170587;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.1pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cs

Note: This product is for in vitro research use only