

Mouse Csnk2b/Casein kinase II subunit beta ELISA Kit

Catalog No: #EK19790

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Csnk2b/Casein kinase II subunit beta ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Casein kinase II subunit beta
Target Name	Casein kinase II subunit beta
Other Names	CK II beta, Phosvitin, Ck2n
Accession No.	P67871
Uniprot	P67871
GeneID	13001;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.34ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Csnk2b

Note: This product is for in vitro research use only