

Rat Csnk2b/Casein kinase II subunit beta ELISA Kit

Catalog No: #EK19791

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Csnk2b/Casein kinase II subunit beta ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Casein kinase II subunit beta
Target Name	Casein kinase II subunit beta
Other Names	CK II beta, Phosvitin, Ck2n
Accession No.	P67874
Uniprot	P67874
GeneID	81650;

Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.32ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Csnk2b

Note: This product is for in vitro research use only