

# Mouse Angptl3/Angiopoietin-related protein 3 ELISA Kit

Catalog No: #EK19792

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Angptl3/Angiopoietin-related protein 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiopoietin-related protein 3
Target Name	Angiopoietin-related protein 3
Other Names	Angiopoietin-like protein 3
Accession No.	Q9R182
Uniprot	Q9R182
GeneID	30924;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.414ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Angptl3

Note: This product is for in vitro research use only