

## KCNN1 Antibody

Catalog No: #44687



Package Size: #44687-1 50ul #44687-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	KCNN1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	KCNN1 antibody detects endogenous levels of total KCNN1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human KCNN1
Target Name	KCNN1
Other Names	hSK1 antibody KCa2.1 antibody KCNN 1 antibody Kcnn1 antibody KCNN1_HUMAN antibody Potassium intermediate/small conductance calcium activated channel subfamily N member 1 antibody SK antibody SK1 antibody SKCa 1 antibody SKCA1 antibody Small conductance calcium activated potassium channel protein 1 a
Accession No.	Swiss-Prot#: Q92952 NCBI Gene ID:3780
Uniprot	Q92952
GeneID	3780;
Calculated MW	60 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only