

## Human CPA3/Mast cell carboxypeptidase A ELISA Kit

Catalog No: #EK19793

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Human CPA3/Mast cell carboxypeptidase A ELISA Kit          |
| Brief Description  | ELISA Kits   |
| Applications       | ELISA  |
| Species Reactivity | Hu   |
| Specificity        | Natural and recombinant Human Mast cell carboxypeptidase A |
| Target Name        | Mast cell carboxypeptidase A                               |
| Other Names        | MC-CPA, Carboxypeptidase A3                                |
| Accession No.      | P15088   |
| Uniprot            | P15088   |
| GeneID             | 1359;  |

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CPA3

Note: This product is for in vitro research use only