

Mouse Cpa3/Mast cell carboxypeptidase A ELISA Kit

Catalog No: #EK19794

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cpa3/Mast cell carboxypeptidase A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Mast cell carboxypeptidase A
Target Name	Mast cell carboxypeptidase A
Other Names	MC-CPA, Carboxypeptidase A3
Accession No.	P15089
Uniprot	P15089
GeneID	12873;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cpa3

Note: This product is for in vitro research use only