

Rat Cpa3/Mast cell carboxypeptidase A ELISA Kit

Catalog No: #EK19795

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cpa3/Mast cell carboxypeptidase A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Mast cell carboxypeptidase A
Target Name	Mast cell carboxypeptidase A
Other Names	Carboxypeptidase A3,R-CPA,RMC-CP
Accession No.	P21961
Uniprot	P21961
GeneID	54242;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cpa3

Note: This product is for in vitro research use only