

Mouse Cd163/Scavenger receptor cysteine-rich type 1 protein M130 ELISA Kit

Catalog No: #EK19796

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd163/Scavenger receptor cysteine-rich type 1 protein M130 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Scavenger receptor cysteine-rich type 1 protein M130
Target Name	Scavenger receptor cysteine-rich type 1 protein M130
Other Names	CD163,M130
Accession No.	Q2VLH6
Uniprot	Q2VLH6
GeneID	93671;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.162ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd163

Note: This product is for in vitro research use only