

Mouse Ctsg/Cathepsin G ELISA Kit

Catalog No: #EK19797

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ctsg/Cathepsin G ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cathepsin G
Target Name	Cathepsin G
Other Names	Vimentin-specific protease,VSP
Accession No.	P28293
Uniprot	P28293
GeneID	13035;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.165ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ctsg

Note: This product is for in vitro research use only