

Human CASP4/Caspase-4 ELISA Kit

Catalog No: #EK19799

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CASP4/Caspase-4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Caspase-4
Target Name	Caspase-4
Other Names	CASP-4,ICE and Ced-3 homolog 2,ICH-2,ICE(rel)-II,Mih1,Protease TX,ICH2
Accession No.	P49662
Uniprot	P49662
GeneID	837;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.08ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CASP4

Note: This product is for in vitro research use only