

Mouse Casp4/Caspase-4 ELISA Kit

Catalog No: #EK19800

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Casp4/Caspase-4 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Caspase-4 |
| Target Name | Caspase-4 |
| Other Names | CASP-4,Caspase-11,CASP-11,Protease ICH-3,Casp11,Caspl,lch3 |
| Accession No. | P70343 |
| Uniprot | P70343 |
| GeneID | 12363; |

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.05ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Casp4

Note: This product is for in vitro research use only