

# Human CDC23/Cell division cycle protein 23 homolog ELISA Kit



Catalog No: #EK19801

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CDC23/Cell division cycle protein 23 homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cell division cycle protein 23 homolog
Target Name	Cell division cycle protein 23 homolog
Other Names	Anaphase-promoting complex subunit 8,APC8,Cyclosome subunit 8,ANAPC8
Accession No.	Q9UJX2
Uniprot	Q9UJX2
GeneID	8697;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.083ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CDC23

Note: This product is for in vitro research use only