

Mouse Cd7/T-cell antigen CD7 ELISA Kit

Catalog No: #EK19802

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cd7/T-cell antigen CD7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-cell antigen CD7
Target Name	T-cell antigen CD7
Other Names	CD7
Accession No.	P50283
Uniprot	P50283
GeneID	12516;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000ng/mL Sensitivity: 3.9ng/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd7

Note: This product is for in vitro research use only