

# Human COLGALT1/Procollagen galactosyltransferase 1 ELISA Kit



Catalog No: #EK19803

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human COLGALT1/Procollagen galactosyltransferase 1 ELISA Kit  |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Hu  |
| Specificity        | Natural and recombinant Human Procollagen galactosyltransferase 1   |
| Target Name        | Procollagen galactosyltransferase 1   |
| Other Names        | Collagen beta(1-O)galactosyltransferase 1,Glycosyltransferase 25 family member 1,Hydroxylysine galactosyltransferase 1,GLT25D1,PSEC0241 |
| Accession No.      | Q8NBJ5  |
| Uniprot            | Q8NBJ5  |
| GeneID             | 79709;  |

## Application Details

Assay Type: Sandwich Detect Range:31.2-2000pg/mL Sensitivity:15.62pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:COLGALT1

Note: This product is for in vitro research use only