

## Human C2/Complement C2 ELISA Kit

Catalog No: #EK19805

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human C2/Complement C2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Complement C2
Target Name	Complement C2
Other Names	C3/C5 convertase
Accession No.	P06681
Uniprot	P06681
GeneID	717;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.086ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:C2

Note: This product is for in vitro research use only