

## Mouse Cd5/T-cell surface glycoprotein CD5 ELISA Kit

Catalog No: #EK19807

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cd5/T-cell surface glycoprotein CD5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse T-cell surface glycoprotein CD5
Target Name	T-cell surface glycoprotein CD5
Other Names	Lymphocyte antigen 1, Ly-1, CD5, Ly-1
Accession No.	P13379
Uniprot	P13379
GeneID	12507;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.421ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cd5

Note: This product is for in vitro research use only