

## Rat Cd5/T-cell surface glycoprotein CD5 ELISA Kit

Catalog No: #EK19808

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cd5/T-cell surface glycoprotein CD5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat T-cell surface glycoprotein CD5
Target Name	T-cell surface glycoprotein CD5
Other Names	CD5
Accession No.	P51882
Uniprot	P51882
GeneID	54236;

## Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.335ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cd5

Note: This product is for in vitro research use only