

Bovine CA1/Carbonic anhydrase 1 ELISA Kit

Catalog No: #EK19809

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CA1/Carbonic anhydrase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Carbonic anhydrase 1
Target Name	Carbonic anhydrase 1
Other Names	Carbonate dehydratase I, Carbonic anhydrase I, CA-I
Accession No.	Q1LZA1
Uniprot	Q1LZA1
GeneID	510801;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20 ng/mL Sensitivity: 0.174 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: CA1

Note: This product is for in vitro research use only