

# Mouse Adamts4/A disintegrin and metalloproteinase with thrombospondin motifs 4 ELISA Kit



Catalog No: #EK19811

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Adamts4/A disintegrin and metalloproteinase with thrombospondin motifs 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse A disintegrin and metalloproteinase with thrombospondin motifs 4
Target Name	A disintegrin and metalloproteinase with thrombospondin motifs 4
Other Names	ADAM-TS 4,Aggrecanase-1
Accession No.	Q8BNJ2
Uniprot	Q8BNJ2
GeneID	240913;

## Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.093ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Adamts4

Note: This product is for in vitro research use only