

## Mouse Bace2/Beta-secretase 2 ELISA Kit

Catalog No: #EK19812

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bace2/Beta-secretase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta-secretase 2
Target Name	Beta-secretase 2
Other Names	Aspartyl protease 1,ASP1,Beta-site amyloid precursor protein cleaving enzyme 2,Beta-site APP cleaving enzyme 2,Memapsin-1,Membrane-associated aspartic protease 1,Theta-secretase
Accession No.	Q9JL18
Uniprot	Q9JL18
GeneID	56175;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:40.7pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bace2

Note: This product is for in vitro research use only