

Mouse Col18a1/Collagen alpha-1(XVIII) chain ELISA Kit



Catalog No: #EK19813

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Col18a1/Collagen alpha-1(XVIII) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(XVIII) chain
Target Name	Collagen alpha-1(XVIII) chain
Accession No.	P39061
Uniprot	P39061
GeneID	12822;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 17pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Col18a1

Note: This product is for in vitro research use only