

# Mouse Adam9/Disintegrin and metalloproteinase domain-containing protein 9 ELISA Kit



Catalog No: #EK19814

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Adam9/Disintegrin and metalloproteinase domain-containing protein 9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 9
Target Name	Disintegrin and metalloproteinase domain-containing protein 9
Other Names	ADAM 9, Meltrin-gamma, Metalloprotease/disintegrin/cysteine-rich protein 9, Myeloma cell metalloproteinase, Kiaa0021, Mdc9, Mltng
Accession No.	Q61072
Uniprot	Q61072
GeneID	11502;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.23ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Adam9

Note: This product is for in vitro research use only