

Mouse Cfd/Complement factor D ELISA Kit

Catalog No: #EK19815

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cfd/Complement factor D ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement factor D
Target Name	Complement factor D
Other Names	28 kDa adipocyte protein, Adipsin, C3 convertase activator, Properdin factor D, Adn, Df
Accession No.	P03953
Uniprot	P03953
GeneID	11537;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.18ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cfd

Note: This product is for in vitro research use only