

Pig CFD/Complement factor D ELISA Kit

Catalog No: #EK19816

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig CFD/Complement factor D ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Complement factor D
Target Name	Complement factor D
Other Names	Adipsin,C3 convertase activator,Properdin factor D,DF
Accession No.	P51779
Uniprot	P51779
GeneID	396877;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/ml Sensitivity:40.9pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CFD

Note: This product is for in vitro research use only