

Human CHR3/Muscarinic acetylcholine receptor M3 ELISA Kit



Catalog No: #EK19817

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CHR3/Muscarinic acetylcholine receptor M3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Muscarinic acetylcholine receptor M3
Target Name	Muscarinic acetylcholine receptor M3
Accession No.	P20309
Uniprot	P20309
GeneID	1131;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 41pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CHR3

Note: This product is for in vitro research use only