



Rat Chrm3/Muscarinic acetylcholine receptor M3 ELISA Kit

Catalog No: #EK19818

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Chrm3/Muscarinic acetylcholine receptor M3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Muscarinic acetylcholine receptor M3
Target Name	Muscarinic acetylcholine receptor M3
Other Names	Chrm-3
Accession No.	P08483
Uniprot	P08483
GeneID	24260;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:41.1pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Chrm3

Note: This product is for in vitro research use only