

Human CD209/CD209 antigen ELISA Kit

Catalog No: #EK19819

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CD209/CD209 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human CD209 antigen
Target Name	CD209 antigen
Other Names	C-type lectin domain family 4 member L,Dendritic cell-specific ICAM-3-grabbing non-integrin 1,DC-SIGN,CD209,CLEC4L
Accession No.	Q9NNX6
Uniprot	Q9NNX6
GeneID	30835;

Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.34ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CD209

Note: This product is for in vitro research use only