

Mouse Bsg/Basigin ELISA Kit

Catalog No: #EK19820

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bsg/Basigin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Basigin
Target Name	Basigin
Other Names	Basic immunoglobulin superfamily,HT7 antigen,Membrane glycoprotein gp42,CD147
Accession No.	P18572
Uniprot	P18572
GeneID	12215;

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.88ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bsg

Note: This product is for in vitro research use only