

Rat Bsg/Basigin ELISA Kit

Catalog No: #EK19821

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bsg/Basigin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Basigin
Target Name	Basigin
Other Names	Glycoprotein CE9,OX-47 antigen,CD147
Accession No.	P26453
Uniprot	P26453
GeneID	25246;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.29ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bsg

Note: This product is for in vitro research use only