

Bovine CSTB/Cystatin-B ELISA Kit

Catalog No: #EK19822

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CSTB/Cystatin-B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cystatin-B
Target Name	Cystatin-B
Other Names	Stefin-B,CST6
Accession No.	P25417
Uniprot	P25417
GeneID	512805;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.064ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CSTB

Note: This product is for in vitro research use only