

Human CSTB/Cystatin-B ELISA Kit

Catalog No: #EK19823

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CSTB/Cystatin-B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cystatin-B
Target Name	Cystatin-B
Other Names	CPI-B,Liver thiol proteinase inhibitor,Stefin-B,CST6,STFB
Accession No.	P04080
Uniprot	P04080
GeneID	1476;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CSTB

Note: This product is for in vitro research use only